

# DRYING OVEN OVE12-225T



# DRYING OVEN OVE12-225T

Heating Oven, Laboratory Drying Oven

Used in Drying, Baking, Disinfecting, Sterilization, Laboratory, Research, Engineering, Industry, Thermostat storage, Heating Storage.

Also known as Heating Oven, Laboratory Drying Oven.

## OVE12-225T DRYING OVEN

Digital display all 304 stainless steel inner chamber and outer

Intelligent PID control system

Large observation window

Tight buckle handle



## SPECIFICATIONS

Model	OVE12-225T
Capacity	225 L
Temperature Accuracy	$\pm 1^{\circ}\text{C}$
Temperature Fluctuation	$\pm 1\%$
Timer Range	0-999 min
Internal Dimension	600Wx500Dx750H mm
Exterior Dimension	965Wx770Dx980H mm
Weight	94/103 kg
Power	2500 W
Power Supply	220 V, 50 Hz

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)