

# THERMOSTATIC DRYING OVEN OVE13-130



# THERMOSTATIC DRYING OVEN OVE13-130

Heating Oven, Laboratory Drying Oven

Used in Drying, Baking, Disinfecting, Sterilization, Laboratory, Research, Engineering, Industry, Thermostat storage, Heating Storage.

Also known as Heating Oven, Laboratory Drying Oven.

## OVE13-130 THERMOSTATIC DRYING OVEN

ALLSENSTM programmable PID control

ALLHEATTM cavity warm-up technology

ALLFLOWTM perfect air current cycling

Pro-INSULATIONTM Insulation technology

SUS304 stainless steel internal material

Build-in memory menu, real-time monitoring

Display parameters: Temperature, Air velocity, Run time, Run/Stop



## SPECIFICATIONS

|                         |                          |
|-------------------------|--------------------------|
| Model                   | OVE13-130                |
| Capacity                | 130 L                    |
| Temperature Range       | RT+5 - 250°C             |
| Ambient Temperature     | 10-30°C                  |
| Temperature Accuracy    | 0.1°C                    |
| Temperature Fluctuation | ±1°C(50 - 240°C)         |
| Temperature Uniformity  | 2%(50 - 240°C)           |
| Ambient Humidity        | <70%                     |
| Timer Range             | 0-99 h or 0-9999 min     |
| Convection Mode         | Forced Convection        |
| Internal Dimension      | 550Wx550Dx450H mm        |
| Exterior Dimension      | 770Wx780Dx693H mm        |
| Package Size            | 890x860x780 mm           |
| Shelves                 | 2                        |
| Weight                  | 65/69 kg                 |
| Power                   | 1270 W                   |
| Power Supply            | Single phase AC220V/50Hz |



**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@m.labstac.com](mailto:contact@m.labstac.com) | Website: [labstac.com](http://labstac.com)