

# VACUUM DRYING OVEN OVE21-25



# VACUUM DRYING OVEN OVE21-25

Vacuum ovens are very versatile pieces of equipment with applications in laboratory research, engineering, and industry. The low pressure environment also minimizes oxidation during drying. A rugged, high quality pressure chamber, gasket seals, and pump are features to look for, as well as convenient, programmable controls and a computer interface.

Used in Biochemical, Pharmaceutical, Healthcare, Agricultural, Environmental protection, disinfection, Sterilization, Drying, Roasting, Laboratory, Research, engineering, Industry,.

Also known as Vacuum Oven, Laboratory Vacuum Oven.

## OVE21-25 VACUUM DRYING OVEN

Stainless steel inner chamber without scroll rolled

Silicone rubber door seal

Strong sealing

Nice design



## SPECIFICATIONS

Model	OVE21-25
Capacity	25 L
Temperature Accuracy	$\pm 1^{\circ}\text{C}$
Temperature Fluctuation	$\pm 2\%$
Vacuum Degree	$< 133 \text{ Pa}$
Timer Range	0-999 min
Internal Dimension	300Wx300Dx280H mm
Shelves/Trays	2
Power	900 W
Power Supply	220 V, 50 Hz

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)