

# TIPS BOX FOR 5000?L TIPS PIC64-05000



# TIPS BOX FOR 5000µL TIPS PIC64-05000

It offers quick and easy access to pipette tips, saving valuable time. Space-saving and environmentally compatible. Commonly used in liquid handling and molecular biology.

Used in Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

## PIC64-05000 TIPS BOX FOR 5000µL TIPS

Made from medical grade Polypropylene

Slippy surface with excellent design

DNase-free

Rnase-free

Non-pyrogenic

Filter and non-filter for choice

Bag and rack are available to pack



## SPECIFICATIONS

Model	PIC64-05000
Volume	5000 µl
Unit	20

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)