

SEROLOGICAL PIPETTE PIC71-1400C



SEROLOGICAL PIPETTE PIC71-1400C

Manufactured for accurate and rapid fluid dispensing. It is sterile, pyrogen-free, cotton plugged, and can be used for a variety of liquid handling applications. The reservoir eliminates the need for multiple reagent containers and saves bench top space.

Used in Plasmid Purification, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

Also known as Laboratory Serological Pipette.

PIC71-1400C SEROLOGICAL PIPETTE

Easy to use

Medical grade original Polystyrene


RNase- free, DNase- free and Non-pyrogenic

Seven types of volume 1ml, 2ml, 5ml, 10ml, 25ml, 50ml and 100ml

Bidirectional graduation improve accuracy

Seven types of color code on tips make pipette

Sterilized by Gamma irradiation



SPECIFICATIONS

Model	PIC71-1400C
Volume	50 ml
Precision	2 %
Negative Graduation (ml)	-3
Graduation	1
Color	Black
Sterile	Gamma
Package	Individual
Pcs/bag	25 pcs/bag
Pcs/ctn	1400 pcs/carton
Unit	1

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com