

# MOTORIZED PIPETTE CONTROLLER PIP32-100B



# MOTORIZED PIPETTE CONTROLLER PIP32-100B

Pipette Controller offers an unmatched control of pipetting speeds from drop by drop or gravity dispense to fast liquid displacement. Designed for comfortable and uncomplicated pipetting even during long lasting sessions. It is light weight and easy to use in laboratory and for research purpose.

Used in Liquid Handling Task, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research.

Also known as Pipetaid, Pipet Gun, Laboratory Pipette Controller.

## PIP32-100B MOTORIZED PIPETTE CONTROLLER

Chemical resistance

Human engineering design

Lightweight and Structure strong

Lightweight and dexterous for operation

A wide range volume for liquid controlling



### SPECIFICATIONS

Model	PIP32-100B
Volume	0.1-100 ml
Colour	Dark blue

### OPTIONAL ACCESSORIES

Accessory Code	Name	MOQ
5200815005	RECHARGING STAND WALL MOUNT BRACKET	20
5200815006	CHARGER American plug	20
5200815007	CHARGER Europe plug	20
5200815008	CHARGER Australian plug	20
5200815009	EACH MEMBRANE FILTERS 0.45UM	20
5200815010	Battery	20

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)