

LOW TEMPERATURE PUMPABILITY & GELATION INDEX TESTER PTL21-A




LOW TEMPERATURE PUMPABILITY & GELATION INDEX TESTER PTL21-A

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

PTL21-A LOW TEMPERATURE PUMPABILITY & GELATION INDEX TESTER

Automatic calibration procedure, graphic display. Multi-point temperature calibration, can be added by user. It has tube locking function, can avoid influencing accuracy when tube is moved for high viscosity.



SPECIFICATIONS

Model	PTL21-A
Lifting method	One-click, Automatic
Viscosity and Gelation curve	Automatic scan
Gelation index and Temperature analysis	Automatic
Temperature	Multi-point calibration

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com