

# SHEAR STABILITY TESTER PTL22-A



# SHEAR STABILITY TESTER PTL22-A

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

## PTL22-A SHEAR STABILITY TESTER

Independent desktop controller, high integrated function button easy operation.  
Automatic stroke counter can cut off power automatically.



## SPECIFICATIONS

Model	PTL22-A
Timing method	Automatic stroke timer
Temperature control range	30 – 35 °C
Rotation speed	925 ± 25 r/min
Power	1100 W

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@m.labstac.com](mailto:contact@m.labstac.com) | Website: [labstac.com](http://labstac.com)