

FLASH POINT AND FIRE POINT TESTER PTL25-A



FLASH POINT AND FIRE POINT TESTER PTL25-A

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

PTL25-A FLASH POINT AND FIRE POINT TESTER

Use heat resistant glass interlayer on the heater,avoid oil splash and fire breaking out.
Automatic stop if there isn't operation after heating 20 minutes. Automatic air cooling
when temperature is too high.



SPECIFICATIONS

| | |
|-------------|--|
| Model | PTL25-A |
| Air cooling | Automatic when temperature is too high |
| Stop | Automatic if there isn't operation after heating 20 minutes |
| Heating | Can only be started if press switch 3 seconds,avoiding touch by mistake. |
| Safety | Use heat resistant glass interlayer on the heater,avoid oil splash and fire breaking out |

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com