

RUST PREVENTION TESTER PTL28-A



RUST PREVENTION TESTER PTL28-A

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

PTL28-A RUST PREVENTION TESTER

Independent quiet stirrer, can test one or several samples. Each stirrer controls rotating speed, test time and automatically stop when test is finished.



SPECIFICATIONS

Model	PTL28-A
Working units	4
Stirrer	Independent quiet, controls rotating speed, test time and automatically stop
Liquid circulator	Sputtering multi-channel
Mixing blade	Sliding automatic
Safety protection device	Prevent over-low liquid level

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@m.labstac.com | Website: labstac.com