

THERMAL STABILITY TESTER PTL29-A



THERMAL STABILITY TESTER PTL29-A

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

PTL29-A THERMAL STABILITY TESTER

Eco-friendly metal bath heating,high effect of heat insulation. Air filling and releasing device with pressure gauge,adjusted by needle valve.



SPECIFICATIONS

Model	PTL29-A
Working units	1 – 4
Temperature	Upto 400 °C
Metal bath heating	Eco-friendly,high effect of heat insulation
Test bomb	Standard stainless steel with inlet valve of 20m1
Safety protection device	Air filling and releasing device with pressure gauge,adjusted by needle valve.

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@m.labstac.com | Website: labstac.com