

AIR RELEASE PROPERTIES TESTER PTL2A-01



AIR RELEASE PROPERTIES TESTER PTL2A-01

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

PTL2A-01 AIR RELEASE PROPERTIES TESTER

Temperature of bath and air are controlled separately. Automatically control blowing time and defoaming time, accuracy can reach 0.01s. Press-on joint is easy to dismantle tube.



SPECIFICATIONS

Model	PTL2A-01
Accuracy	Upto 0.01 s
Samples size	2, simultaneously
Tube	Press-on joint make it easy to dismantle
Blowing time and Defoaming time control	Automatic
Liquid circulator	Sputtering multi-channel
Temperature of bath and air	Controlled separately
Screen	Display pressure
Safety	Safety protection device an prevent over-low liquid level, with needle valve to adjust

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com