

COPPER STRIP CORROSION TESTER PTL2B-06



COPPER STRIP CORROSION TESTER PTL2B-06

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

PTL2B-06 COPPER STRIP CORROSION TESTER

Can test lubricating oil and grease with different cylinders. Eco-friendly metal bath heating, no need heat-conducting medium. There is prompt when test is finished.



SPECIFICATIONS

Model	PTL2B-06
Samples size	6, simultaneously
Test samples	Lubricating oil and grease with different cylinders
Safety	Safety protection device an prevent over-high liquid level
Metal bath	Eco-friendly heating, no need heat-conducting medium

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com