

# WATER SEPARABILITY TESTER PTL2C-09



# WATER SEPARABILITY TESTER PTL2C-09

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

## PTL2C-09 WATER SEPARABILITY TESTER

Sliding automatic location device, no need to adjust blade. Sputtering scroll liquid circulator can guarantee temperature uniformity of each point very well.



### SPECIFICATIONS

Model	PTL2C-09
Working units	9
Stirrer	Digital, rotating, adjustable speed
Stirring time and record separation time	Automatically controlled
Unit of separation time	Selective
Safety	Safety protection device an prevent over-low liquid level
Liquid circulator	Sputtering scroll

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)