

DENSITY TESTER PTL2D-02



DENSITY TESTER PTL2D-02

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

PTL2D-02 DENSITY TESTER

Low noise magnetic pump cross circulation structure. Can test density of antifreeze



SPECIFICATIONS

Model	PTL2D-02
Working units	2
Test sample	Antifreeze
Temperature control range	Ambient to 20 °C ± 0.1 °C
Refrigeration	Shock absorption compressor, fluorine free
Bath	Low temperature with lighting
Magnetic pump	Low noise with cross circulation structure
Safety	Safety protection device to prevent over-low liquid level

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com