

FREEZING POINT TESTER PTL2K-02




FREEZING POINT TESTER PTL2K-02

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

PTL2K-02 FREEZING POINT TESTER

Make two test simultaneously, calculate average value. Replace magnetic stirrer with mechanical stirrer, which can avoid uneven mixing and influencing measuring accuracy. Adjustable tube clamp can fix test tube perfectly.



Product Image Coming Soon

SPECIFICATIONS

Model	PTL2K-02
Samples size	2 sample simultaneously
Calculate	Average value
Temperature	- 70 °C
Refrigerant	Dual air compressors
Stirrer	Mechanical stirrer, which can avoid uneven mixing and influencing measuring accuracy
tube clamp	Adjustable can fix test tube perfectly
Lowered temperature curve	Automatic
Freezing point calculation	Automatic

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com