

POUR POINT TESTER PTL2L-02



POUR POINT TESTER PTL2L-02

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

PTL2L-02 POUR POINT TESTER

Dry type over speed refrigeration technology without medium. Temperature of two troughs are controlled separately. Compact structure, can put on platform directly



SPECIFICATIONS

| | |
|-------------------|---|
| Model | PTL2L-02 |
| Trough | Two, temperature controlled separately |
| Samples size | 2 samples in each trough |
| Temperature Range | Reach -30°C in 30 minutes and the lowest temperature is -70°C , |
| Refrigeration | Dry type over speed refrigeration technology without medium |
| Refrigerant | Fluorine free |
| Structure | Compact structure, can put on platform directly |

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com