

BOILING POINT TESTER PTL2M-02



BOILING POINT TESTER PTL2M-02

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

PTL2M-02 BOILING POINT TESTER

Comply with test requirements of antifreeze and brake fluid. Digitally display control voltage. High temperature resistant glass partition



SPECIFICATIONS

Model	PTL2M-02
Furnace	IR heating furnace no radiation, fast temperature rising speed
Test requirements	Comply with test requirements of antifreeze and brake fluid
Glass partition	High temperature resistant
Display	Control voltage

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com