

CORROSION- PREVENTIVE PROPERTIES TESTER PTL2N-A



CORROSION-PREVENTIVE PROPERTIES TESTER

PTL2N-A

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

PTL2N-A CORROSION-PREVENTIVE PROPERTIES TESTER

Bearing seat made of engineering plastics can avoid corrosive influence caused by metal bearing seat. Automatically stop motor when test is finished. Can do static and dynamic tests.



SPECIFICATIONS

Model	PTL2N-A
Tests	Static and Dynamic
Motor	Automatic stop when test is finished
Operation	Observe bearing through reflector with zoom function
Time	Can edit working time and stop time
Bearing seat	Made of engineering plastics can avoid corrosive influence

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com