

# CORROSION TESTER PTL20-03



# CORROSION TESTER PTL2O-03

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

## PTL2O-03 CORROSION TESTER

Sputtering multi-channel liquid circulator can guarantee temperature uniformity very well. Condenser quick detachable device



## SPECIFICATIONS

Model	PTL2O-03
Working units	3, can be extended to 4-6 units
Circulator	Sputtering multi-channel liquid circulator can guarantee temperature uniformity
Condenser	Quick detachable device
Safety protection device	Can prevent over-low liquid level

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)