

WATER WASHOUT RESISTANCE TESTER PTL2P-60



WATER WASHOUT RESISTANCE TESTER PTL2P-60

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

PTL2P-60 WATER WASHOUT RESISTANCE TESTER

Unique design, no need to start test shaft, and start water pump circulating system independently. Upgraded connecting shaft can quickly install and dismantle bearing and inside and outside baffle.



SPECIFICATIONS

| | |
|-------------------|---|
| Model | PTL2P-60 |
| Motor | fixed speed 600 rpm |
| Temperature | Control test temperature at 38 °C and 79 °C as per ASTM standard |
| Independent start | Test shaft and Start water pump circulating system |
| Shaft | Upgraded connecting shaft can quickly install and dismantle bearing and inside and outside baffle |

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com