

OXIDATION STABILITY TESTER PTL2R-04



OXIDATION STABILITY TESTER PTL2R-04

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

PTL2R-04 OXIDATION STABILITY TESTER

Can re-evaluate every test curve. By setting the evaluating time range, take certain section of curve to evaluate, which can avoid some sections with abnormal data caused by washing, sample or other reasons. Independent flow control system can automatically control flow of every sample.



SPECIFICATIONS

| | |
|--|---|
| Model | PTL2R-04 |
| Samples size | 4 sample simultaneously |
| Flow control system | Independent system, can automatically control flow of every sample. |
| Curve of conductivity and Induced knee point | Automatically recorded |
| Test curve | Can re-evaluate every test curve |
| Software | Functions of multi-curve comparing and curve analysis, measure induced response time and evaluate stable time |

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com