

COLD FILTER PLUGGING POINT TESTER PTL2S-05



COLD FILTER PLUGGING POINT TESTER PTL2S-05

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

PTL2S-05 COLD FILTER PLUGGING POINT TESTER

Only need to set expected temperature of cold filter plugging point, button start and finish test automatically. Automatically release pressure in tube after 1 min suction filtration

SPECIFICATIONS

Model	PTL2S-05
Detection	Non-contact IR detection technology
Vacuum tanks	Two 5L
Temperature	Only need to set expected temperature of cold filter
Plugging point, Button start and Finish test	Automatic
Pressure	Automatic pressure release in tube after 1 min suction filtration

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com