

# VAPOUR PRESSURE TESTER PTL2T-02



# VAPOUR PRESSURE TESTER PTL2T-02

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

## PTL2T-02 VAPOUR PRESSURE TESTER

Equipped with quick connector, make sure the test can be performed quickly and easily. Automatic temperature control, alarming while reaching required test temperature.

### SPECIFICATIONS

Model	PTL2T-02
Test units	Two test units allow independent test for two samples at the same time
Display	Big colorful touch screen is convenient to use
Test accuracy	360° horizontal flip and oscillation improves test accuracy
Temperature control	Automatic temperature control, alarming while reaching required test temperature
Connector	Equipped with quick connector, make sure the test can be performed quickly and easily

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)