

# WATER CONTENT TESTER PTL2Z-15



# WATER CONTENT TESTER PTL2Z-15

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

## PTL2Z-15 WATER CONTENT TESTER

The instrument is suitable to determine the water content of petroleum products and to determine the water content of lubricating grease. It adopts double-unit structure. It can do determination for two samples. The gripper is designed reasonably. The installation and dismantlement are convenient.



## SPECIFICATIONS

Model	PTL2Z-15
Heating control	Can be continuously adjusted by a silicon knob, manual.
Ambient temperature	(15-35) °C
Relative humidity	≤ 85%
Power supply	AC(220±10%)V, 50Hz
Total power consumption	≤2200 W
Dimension	430×320×700 mm
Net weight	8.5kg

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)