

ANILINE POINT TESTER PTL2Z26-02



ANILINE POINT TESTER PTL2Z26-02

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

PTL2Z26-02 ANILINE POINT TESTER

The instrument can meet requirements for determining aniline point of transparent oils and aniline point of dark color oils. Digital temperature controller. High temperature control precision. Stable and reliable. Easy to operate.



SPECIFICATIONS

Model	PTL2Z26-02
Temperature range	Ambient to 200°C
Temperature precision	≤0.5°C
Heating power	≤2100W
Ambient temperature	15°C~35°C
Relative humidity	≤ 85%
Power supply	AC(220±10%)V, 50Hz
Dimension	370×520×460 mm
Net weight	14 kg

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com