

MODERATELY HIGH TEMPERATURE PISTON DEPOSITS APPARATUS PTL41-05



MODERATELY HIGH TEMPERATURE PISTON DEPOSITS APPARATUS PTL41-05

The instrument is used to test thermo-oxidation engine oil simulation. The procedure includes determination of mass of deposit formed on a specially constructed test rod after oxidative and catalytic conditions at 285 °C is repetitively exposed to passage of 8.5 g of engine oil.

Used in Petroleum Industry, Oil and Gas Industry.

PTL41-05 MODERATELY HIGH TEMPERATURE PISTON DEPOSITS APPARATUS

User can set working time, automatic stop after finishing test. Catalyst magnetic mixing device can set mixing time and speed. All parameters can be set on touch screen, one-click start test.



Product Image Coming Soon

SPECIFICATIONS

Model	PTL41-05
Gas flow	Electronic flow controller accurately controls the gas flow at 10ml±0.2ml/min
Mixing device	Catalyst magnetic mixing device can set mixing time and speed
Display	All parameters can be set on touch screen, one-click start test
Optional	ASTM D6335 and SH/T0750 test components

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@m.labstac.com | Website: labstac.com