

# PORTABLE REFRIGERATOR REF31-240



# PORTABLE REFRIGERATOR REF31-240

A portable refrigerator is self-contained storage chest which uses an outside power source to control its internal temperature, also it doesn't rely on ice or freezer packs to keep its contents cool.

Used in Industry, Hospitals ,Cosmetic, Pharmaceutical, Electronics, Laboratory, Medical, Research.

## REF31-240 PORTABLE REFRIGERATOR



### SPECIFICATIONS

|                       |  |
|-----------------------|--|
| Model                 | REF31-240  |
| External Size         | 330x235x275 mm   |
| Internal Size         | 240x145x180 mm   |
| Inner Capacity        | 6 L  |
| Product Configuration | 1pc refrigerator, 4 pcs Crystal storage ice box 400ml, 1pc digital thermometer |
| Insulation Material   | PU   |
| Gross Weight          | 2.8 kg   |
| Packing Size          | 380x260x290 mm   |

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)