

INCUBATOR SHAKER SHA024-72B



INCUBATOR SHAKER SHA024-72B

SHA024-72B INCUBATOR SHAKER

Adopt 4.3 inches touch screen, intelligent touch screen control, parameters real time display, operation easily

Temp., rotation speed and time etc parameters fast setting

Combine incubator, shaker functions together with small space, huge capacity

Colorful touch screen controller, menu operation interface, multi datas display at one screen, easily to observe and operate

Fluorine free design

User's set parameters can be stored automatically when power off suddenly, and it resume last program settings when power on

PID microcomputer control temp. and shake frequency with timer

Compressor and cycle fan etc key components are reliable suppliers, environmental non-free refrigerants

Brushless DC motor, large start torque, wide speed adjustment, free maintenance

Sterilization system: The ultraviolet lamp can sterilize the chamber regularly

Protection on instruments: Comply international standard secondary temp. limiter alarm system, alert the operator with sound and light alarms, ensure operator is safe without any accident

Protection on key components: Key components have over current, over temp., over load etc safety protection, it can prevent instruments accidents without precautions

Protection on samples: Working chamber temperature higher or lower the set temp., alarm starter cut down the heater, alert the operator with sound and light alarms

Protection on operator: Cabinet and door is special designed for insulation, low heat of cabinet body, ensure operator use instrument without burnt



SPECIFICATIONS

| | |
|------------------------|-----------------------------|
| Model | SHA024-72B |
| Speed Range | 40-300 rpm |
| Frequency accuracy | ±1 rpm |
| Amplitude | 26 mm |
| Temp. Control Range | 4 - 65°C |
| Temperature Resolution | ±0.1°C |
| Temperature Uniformity | at 37°C±0.8°C |
| Timing Range | 0 - 99h59 min |
| Shelves size | 920x500 mm (double shelves) |
| Overall Dimension | 1200Wx875Dx1370H mm |
| Power | 2250 W |
| Standard configuration | 250mlx90pcs |



Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@labstac.com | Website: labstac.com