

# INCUBATOR SHAKER SHA027-700A



# INCUBATOR SHAKER SHA027-700A

## SHA027-700A INCUBATOR SHAKER

Large LCD screen to display more data at same time

R134a refrigerant, imported compressor and fan motor

Big observation windows

304 Stainless steel chamber and platform, easy to clean.

There is a 25mm instruction connection hole on the left side of the chamber for easy testing operation and temperature measurement

The parameters can be automatically stored in case of power failure, and it will continue run as presetting program after turn on

Microprocessor PID controller for temperature and shaking speed with timing function



## SPECIFICATIONS

Model	SHA027-700A
Speed Range	40-300 rpm
Electrical Requirement	220 V 50 H
Amplitude	26 mm
Temperature Range	RT+5 - 65°C
Display Resolution	0.1°C
Timing Range	0~5999 min
Power Consumption	1900 W
Platform Dimension	920x500 mm
Platform quantity	2 pcs

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)