

INCUBATOR SHAKER SHA27-0300



INCUBATOR SHAKER SHA27-0300

In Incubator Shaker, the constant agitation with uniform amplitude provides the optimum oxygen transfer to the cultures, leading to better and consistent yield. A variety of platforms and different laboratory accessories permits configuration to fulfill your demands.

Used in Cell Cultures, Cell Aeration, Microbiology, Increasing Solubility Rates, Metabolism Studies, Bacterial Cultures, Bacteriology.

Also known as Laboratory Incubator Shaker, Incubator Shaker.

SHA27-0300 INCUBATOR SHAKER

Speed can be set as 60-400 rpm

The electronic control of closed cycle adjustable heating system

Sainless steel interior and Electrostatic sprayed exterior

Call to restore function, without power supply

With strong influence of intermittent rapid heating system of.200 type and 100B type with air pressure regulating valve

Flexible and convenient, wide variety of detachable accessories



SPECIFICATIONS

| | |
|------------------------|---------------------------------------|
| Model | SHA27-0300 |
| Speed Range | 20-300 rpm |
| Standard Capacity | 500 mlx3, 250 mlx3, 100 mlx4, 50 mlx4 |
| Temperature | RT+5°C~60°C |
| Temperature Accuracy | ±0.1°C |
| Temperature Uniformity | ±1°C |
| Shaking Speed Accuracy | ±1 rpm |
| Timing Range | Optional |
| Motion | Orbital |
| Amplitude | 24 mm |
| Tray | 1 |
| Platform Dimension | 4375x405 mm |
| Inner Dimension | 428x456x360 mm |
| Overall Dimension | 640x605x520 mm |
| Display | LCD |
| Weight | 72 kg |
| Power | 365 W |
| Power Supply | 220 V±10% |



Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@labstac.com | Website: labstac.com