

3D GYRATORY ROCKER SHA61-20



3D GYRATORY ROCKER SHA61-20

Rocking Shaker are designed to deliver salient performance, accuracy, and reproducible results for all your application needs. Lightweight and portable for easy transfer, outstanding uniform motion, high quality stain resistant platforms with non-slip rubber mats and low noise are the efficient and reliable features that makes it unique.

Used in Laboratory, Research, Mixing, Incubation.

Also known as Laboratory Rocking Shaker.

SHA61-20 3D GYRATORY ROCKER

Gentle but thorough mixing

Suits for blood collection tubes and centrifuge tubes

Tilt angle and speed have been optimally set

Safe for use in cold rooms and incubators



SPECIFICATIONS

| | |
|---------------|---------------------------|
| Model | SHA61-20 |
| Speed Range | 20 rpm (fixed) |
| Shaking Orbit | 3D / 20 (fixed) |
| Max. Load | 0.8 kgs |
| Voltage | AC 220V/AC 110V, 50/60 Hz |
| Power | 45 W |
| Dimension | W.240 x D.170 x H.150 mm |
| Net Weight | 1.6 kgs |

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com