

# THERMO SHAKER INCUBATOR SHA74-1500



# THERMO SHAKER INCUBATOR SHA74-1500

The unique design principle which combines mixing and shaking with intelligent control makes its operation simpler, safe, and convenient. It is compact, highly-efficient and consists of several features including intelligent programming for easy and accurate operation results. It is primarily used to mix different micro capacity tubes, PCR plates, micro-titer plates and other laboratory equipment.

Used in Incubation, Cultivation, Re-suspending Pellets, Denaturation, PCR preparations, ELISA Assays, Colorimetric Assays, Enzyme, Protein Analysis.

Also known as Thermo Shaker, Lab Shaker, Thermo Shaker Incubator, Laboratory Thermo Shaker Incubator.

## SHA74-1500 THERMO SHAKER INCUBATOR

Optimized 2D mixing designed for simultaneous mixing and vortexing applications

Pre-programmed soft keys and reliable universal holder

Automatic imbalance detection and anti-spill technology

Mixing in a 3-in-1 format: Plates (up to 96-well and 384-well formats), Micro-Tubes, vortexing numerous tube formats

Custom blocks available for different experiment needs



## SPECIFICATIONS

Model	SHA74-1500
Speed Range	200-1500 rpm
Temperature	0-100°C
Temperature Accuracy	<± 0.5°C
Temperature Resolution	0.1°C
Vortex Shaking Speed	200-1500 rpm
Shaking Orbit	3 mm
Time of Heating	≤20 min (from 20°C to 100°C)
Time of Cooling	≤25 min (from ambient to ambient -20°C)
Timing Range	1 min~99 h59 min
Motion	Orbital
Overall Dimension	300x220x170 mm
Weight	8.5 kg
Power	150 W
Power Supply	AC220 V 50/60 Hz or AC120 V 50/60 Hz

## OPTIONAL ACCESSORIES

Accessory Code	Name	Description
1901907005	Block	54x0.5 ml
1901907006	Block	35x1.5 ml

Accessory Code	Name	Description
1901907012	Block	24x5.0ml tubes
1901907013	Block	15x5.0 ml Eppendorf tubes

1901907007	Block	35x2.0 ml
1901907008	Block	15x0.5 ml+20x1.5 ml
1901907009	Block	24xDia. ≤φ12 mm tube
1901907010	Block	103x67x30 Rectangular bath block
1901907011	Block	96 wells microplate

1901907014	Block	96x0.2 ml
1901907015	Block	12x15 ml Falcon tube
1901907016	Block	6x50 ml Falcon tube
1901907017	Block	96-Deep well microplate

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
 Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)