

# ATOMIC ABSORPTION SPECTROPHOTOMETE R SPE41-0110



# ATOMIC ABSORPTION SPECTROPHOTOMETER

## SPE41-0110

Used in Food and Beverage Industry, Water Analysis, Clinical Research, Pharmaceutical, Mining and Geology, Environmental Monitoring, Oil and Petroleum, Forensics..

Also known as Atomic Absorption Spectrometry, Metal Analysis Spectroscopy..

## SPE41-0110 ATOMIC ABSORPTION SPECTROPHOTOMETER

Six lamp flame method. Full titanium combustion head, 50 mm or 100 mm general combustion head and polymer explosion-proof spray chamber. With customisable atomizer efficient glass atomizer.



### SPECIFICATIONS

Model	SPE41-0110
Light Source	≤ 8 lamps automatic turret, automatic alignment
Power Supply	110/220 V (+5% ~ -10%), 60/50 Hz; 5000 VA
Lamp Current	Pulsed power supply
Optical System	Large 1800 /mm grating ruling, full closed optical system
Wavelength Range	185 - 900 nm. Automatically peak find,a key opticaloptimization function
Wavelength Accuracy	± 0.15 nm
Wavelength Repeatability	< 0.10 nm
Spectral Bandwidth	0.1, 0.2, 0.4, 0.7, 1.0, 1.4, 2.0 nm (7 steps with automatic changeover)
Baseline Stability	≤ ± 0.002A/30 min (Static) ≤ ± 0.004A/30min (Dynamic)
Absorbance Range	0 - 4 A
Flame Analytical System	
Detector:	Imported photomultiplier tube
Burner Head:	Full titanium combustion head, 50 mm or 100 mm general combustion head
Atomization Chamber:	Polymer explosion-proof spray chamber
Nebulizer:	Atomizer efficient glass atomizer, can also be customized
Ignition Type:	Microcomputer control, automatic ignition
Gas Control:	Automatic gas control system
Detection Limits(Cu):	0.002 µg/mL
Precision:	RSD ≤ 0.5%



**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)