

UVC STERILIZER STE14-1360



UVC STERILIZER STE14-1360

UVC Sterilizer is cost effective, Chemical free, efficient and environment friendly used to kill bacteria, viruses, mold, fungi and micro-organisms.

Used in Clinics, Disinfection of Microorganisms, Research, Pharmaceuticals.

Also known as Germicidal Lamp.


STE14-1360 UVC STERILIZER

Simple, easy and convenient-Connected to a wireless network for operation via a mobile phone or tablet devices

Compact and lightweight design-Easy to setup and flexible for any space

Effective disinfection-By using UVC germicidal lamps but no chemicals, it can effectively kill 99.9% germs, also no need to clean up after use

Schedule and timer-Easy to use a mobile phone or tablet to set up timer or schedule for disinfection



SPECIFICATIONS

Model	STE14-1360
Wavelength	253.7 nm
Lamp Type	UVC Lamp
Intensity	1360 uW/cm ²
Lamp Holder	2G11 Base
Wi-Fi	IEEE 802.11 b/g/n 2.4GHz
Android 4.1, iOS 9.0 or higher version	OS Language
English, French, Italian, German, Bulgarian, Czech, Danish, Greek, Polish, Portuguese, Russian, Slovak, Spanish, Traditional Chinese, Turkish	Dimension (WxHxD) and Weight(kg)
260x170x475 mm 2.9 Kg	
495x235x165 mm 3.3 Kg	Package Dimension and Weight
144 Watt	Power
	Input

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com