

MAGNETIC STIRRER STR12-003



MAGNETIC STIRRER STR12-003

It is distinctly designed to meet the users application. Small, exquisite, convenient, portable, simple operation and long durability makes it an ideal choice.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Magnetic Stirrer, Digital Magnetic Stirrer.

STR12-003 MAGNETIC STIRRER

Speed range of 0-1500 rpm

Maximum stirring quantity 3L

ABS casing with fire retardant material, resistance of weak acid and alkali



SPECIFICATIONS

| | |
|--------------------|---|
| Model | STR12-003 |
| Capacity | 3 L |
| Speed Range | 0~1500 rpm |
| Temperature | 5-40°C |
| Humidity | 0.8 |
| Motor Input/output | 5W/3W |
| Work Plate Size | 135 mm |
| Stirrer Dimension | 50 mm |
| Platform Material | Plastic |
| Speed Display | Scale |
| Motor Type | DC motor |
| Protection Class | IP42 |
| Dimension | 150x260x80 mm Package Size (WxDxH)-220x295x180 mm |
| Weight | 2 kg |
| Power | 10 W |
| Power Supply | AC100/200V±10%, 50/60Hz |

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@mail.labstac.com | Website: labstac.com