

# MAGNETIC STIRRER STR14-008



# MAGNETIC STIRRER STR14-008

It is distinctly designed to meet the users application. Small, exquisite, convenient, portable, simple operation and long durability makes it an ideal choice.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Magnetic Stirrer, Digital Magnetic Stirrer.

## STR14-008 MAGNETIC STIRRER

Speed range of 15-1500rpm

Automatic reverse rotation for better mixing result

High IP protection class (IP 65)

Set-up plate is made from chemically resistant materials

Slip-proof and stable



## SPECIFICATIONS

Model	STR14-008
Capacity	0.8 L
Speed Range	15-1500 rpm
Temperature	5-40°C
Humidity	Permissible Relative Humidity-0.8
Motor Input/output	5W/3W
Maximum Load Capacity	2 kg
Work Plate Size	90x90 mm
Platform Material	ABS
Protection Class	IP65
Dimension	175x120x15 mm
Weight	0.3 kg
Power Supply	110-240V, 50/60 Hz

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)