

# MAGNETIC STIRRER STR15-020



# MAGNETIC STIRRER STR15-020

It is distinctly designed to meet the users application. Small, exquisite, convenient, portable, simple operation and long durability makes it an ideal choice.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Magnetic Stirrer, Digital Magnetic Stirrer.

## STR15-020 MAGNETIC STIRRER

Speed range of 100-1500 rpm

Max. stirring quantity 3L

LED display set and actual speeds

Microprocessor control technology guarantees speed control accuracy

ABS casing with fire retardant material, resistance of weak acid and alkali



## SPECIFICATIONS

Model	STR15-020
Capacity	20 L
Speed Range	100~1500 rpm
Temperature	5-40°C
Humidity	0.8
Motor Input/output	5W/3W
Work Plate Size	135 mm
Stirrer Dimension	50 mm
Platform Material	Plastic
Speed Display	LED
Motor Type	DC motor
Protection Class	IP42
Dimension	150x260x80 mm Package Size (WxDxH)-220x295x180 mm
Weight	2 kg
Power	15 W
Power Supply	AC100/200V±10%, 50/60Hz

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)