

# MAGNETIC STIRRER STR17-010



# MAGNETIC STIRRER STR17-010

It is distinctly designed to meet the users application. Small, exquisite, convenient, portable, simple operation and long durability makes it an ideal choice.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Magnetic Stirrer, Digital Magnetic Stirrer.

## STR17-010 MAGNETIC STIRRER

Speed range of 0-1500 rpm

Maximum stirring quantity 10 L

Glass ceramic plate, excellent chemical resistance and easy to clean



### SPECIFICATIONS

Model	STR17-010
Capacity	10 L
Speed Range	0-1500 rpm
Temperature	5-40°C

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)