

MAGNETIC HOTPLATE STIRRER STR21-000



MAGNETIC HOTPLATE STIRRER STR21-000

Rugged and compact magnetic hotplate stirrers are ergonomically designed and provide uniform heating surfaces. It delivers reliable and consistent stirring quickly and efficiently. The low-profile design of these magnetic hotplate stirrers makes them easy to store and fit into fume hoods.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Magnetic Hotplate Stirrer, Digital Magnetic Stirrer.

STR21-000 MAGNETIC HOTPLATE STIRRER

Available for heating and stirring standard/non-standard reaction flasks between 50 ml and 20 L

DC brushless motor make stable operation, low noise, long using life and no spark

Metal plate shell, firm and durable, thermostability and anti-corrosion

The control panel designed 30° slopes for convenient operation

Steady working at low speed, powerful working at high speed



SPECIFICATIONS

| | |
|--------------------|------------------------|
| Model | STR21-000 |
| Capacity | 20 L |
| Stirring Positions | 1 |
| Speed Range | 0~1800 rpm |
| Timing Range | continuous |
| Stirrer Dimension | 80 mm |
| Platform Material | Silicone Pad |
| Platform Dimension | 154x163 mm |
| Overall Dimension | 156Wx270Dx70H mm |
| Fuse | 250 V,1 A, ϕ 5x20 |
| Weight | 2.6 kg |
| Power | 30 W |
| Power Supply | AC100-230 V, 50/60 Hz |

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com