

# MULTISTATION MAGNETIC STIRRER STR34-10



# MULTISTATION MAGNETIC STIRRER STR34-10

Multistation Magnetic Stirrer is designed for high-efficiency experiments and has stir bar configurations to choose. Laboratory space can be used effectively. Low speed but stable, high speed but powerful with low noise. Non-slip and anti-corrosion.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Multistation Magnetic Stirrer, Digital Magnetic Stirrer.

## STR34-10 MULTISTATION MAGNETIC STIRRER

Speed range of 0-1100 rpm

High performance magnetic stirrer with 10 stirring positions benefit laboratory space

Individually controlled stirring for consistency with various samples

Brushless DC motor, maintenance free and provides constant speed

Stainless steel plate with silicone film, anti-slip and anti-corrosion



## SPECIFICATIONS

Model	STR34-10
Capacity	400 ml
Stirring Positions	10
Speed Range	0~1100 rpm
Temperature	Permissible Ambient Temperature -5-40°C
Permissible Relative Humidity	0.8
Work Plate Size	180x450 mm
Platform Material	Stainless steel with silicon film
Speed Display	Scale
Motor Type	Brushless DC motor
Dimension (WxDxH)	182x552x65 mm
Weight	7 kg
Power Supply	AC110/220V, 50/60Hz

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)