

# MULTISTATION MAGNETIC STIRRER STR37-10



# MULTISTATION MAGNETIC STIRRER STR37-10

Multistation Magnetic Stirrer is designed for high-efficiency experiments and has stir bar configurations to choose. Laboratory space can be used effectively. Low speed but stable, high speed but powerful with low noise. Non-slip and anti-corrosion.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Multistation Magnetic Stirrer, Digital Magnetic Stirrer.

## STR37-10 MULTISTATION MAGNETIC STIRRER

Heating temperature up to 120°C

High performance multi-position magnetic stirrer with uniform temperature distribution

Individually controlled stirring for consistency with various samples

Individually controlled stirring for consistency with various samples

Electronic speed control motor provides constant speed

DC brushless motor maintenance free and low noise

Stainless steel plate with silicone films, anti-slip and anti-corrosion



## SPECIFICATIONS

Model	STR37-10
Capacity	400 ml
Stirring Positions	10
Speed Range	0~1100 rpm
Temperature	RT-120°C Maximum Temperature-140°C Permissible Ambient Temperature-5-40°C
Permissible Relative Humidity	0.8
Work Plate Size	180x450 mm
Platform Material	Stainless steel with silicon film
Speed Display	Scale
Motor Type	DC motor
Dimension (WxDxH)	182x552x65 mm
Weight	7 kg
Power Supply	AC110/220V, 50/60Hz

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)