

LABSTAC

MAGNETIC STIRRER STR48-1500



MAGNETIC STIRRER STR48-1500

It is distinctly designed to meet the users application. Small, exquisite, convenient, portable, simple operation and long durability makes it an ideal choice.

Used in Cell Culture, Laboratory, Research, Medical.

Also known as Laboratory Magnetic Stirrer, Digital Magnetic Stirrer.

STR48-1500 MAGNETIC STIRRER

Ultra-light and Ultra-thin, easy to store and save laboratory space

Advanced non-wearing magnetic coil technology, quieter and more stable operation, longer service life.

Can use mobile power supply, convenient for outdoor use, not limited by experimental conditions

Counter clockwise and clockwise counter clockwise stirring functions can be switched to make

the stirring more uniform and meet a variety of experimental needs.

The LCD screen displays the parameters in real time, the LED indicator shows the working status, and the instrument running status is more intuitive.

Touch buttons and sliding touch adjustment, durable, novel and stylish.

Dust and water proof, high protection level (IP65).



SPECIFICATIONS

Model	STR48-1500
Capacity	Max. stirring amount (H ₂ O)-1 L
Speed Range	20 ~ 1500 rpm
Time Range	0 ~ 9hours (0 is forever)
Speed Display	LCD
Platform Material	PVC
Platform size	φ110
Clockwise/counter clockwise stirring intermittent time (except single direction stirring)	45s/5s
Dimension	110x180x14.5 mm
Weight	0.18 Kg
Power Supply	DC 5V 2A
Protection Class	IP rating-IP63

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com