

HOMOGENIZER STR54-18000



HOMOGENIZER STR54-18000

Homogenizers are used for fast organization homogeneous, separation and emulsification of samples. It is well-suited for a wide range of processing applications. For Bio- pharmaceutical and research field along with safety measures, this product becomes unique choice for the buyers to purchase according to their needs.

Used in Drug Development, Food & Beverage Processing, Pharmaceutical, Laboratory, Scientific, Emulsions, Research. Also known as Laboratory Homogenizer, Rotor Stator Homogenizer, Laboratory Rotor Stator Homogenizer.

STR54-18000 HOMOGENIZER

Light and easy to operate

Use micro DC motor for grinding trace tissue, grinding microstructure, short time and quick effect

Grinding pestle is light and reusable, suitable for various sterile ways, convenient to be dismantled



SPECIFICATIONS

| | |
|---------------------------------|----------------------|
| Model | STR54-18000 |
| Speed Range | 18000rpm Fixed Speed |
| Working Noise | 55dB |
| Battery non loaded Working time | 2 hours (Idle) |
| Battery Charging Time | > 2 Hours |
| Battery | AA battery x2 |
| Dimension | 143x25x143 mm |
| Weight | 0.18 Kg |
| Voltage | 2.4V |

OPTIONAL ACCESSORIES

| Accessory Code | Name | Handling Capacity |
|----------------|------------------------|------------------------|
| 3300810005 | Plastic Pestle | plastic pestle (10pcs) |
| 3300810006 | Stainless Steel Pestle | Stainless Steel Pestle |

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@m.labstac.com | Website: labstac.com