

# DIGITAL DESKTOP ULTRASONIC CLEANER ULC11-1400



# DIGITAL DESKTOP ULTRASONIC CLEANER

## ULC11-1400

Digital ultrasonic cleaner is effective way to clean lab tools. It is a high efficiency cleaner specifically developed to clean surgical instruments.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## ULC11-1400 DIGITAL DESKTOP ULTRASONIC CLEANER

High strength and environmental ABS injection molding

Top quality ultrasonic ceramics parts, quality higher than competitors

Industrial IC, high voltage isolation design, safe using and long working life

Damp proof and anti corrosion process of the PCB

Strict production and quality test control

Digital control, multi stage ultrasonic cleaning adjustment

The tank is polished by SUS304 stainless steel stamping point

2500V high voltage test to ensure the product safety, strict production process, aging process and quality control



## SPECIFICATIONS

Model	ULC11-1400
Capacity	1400 ml
Working frequency	40 Khz
Ultrasonic Power	60 W
Heating power	No
Timing time	1-5 mins
Temperature range	No
Cover	Yes
Drain valve	No
Washing basket	Standard
Inner tank dimension	168x136x62 mm
Packing size	316x313x213 mm
Gross weight	2.25 Kg
Power Supply	220V+ 10% 50Hz/60Hz 110V+ 10% 50Hz/60Hz

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)