

DIGITAL DESKTOP ULTRASONIC CLEANER ULC11-2500



DIGITAL DESKTOP ULTRASONIC CLEANER

ULC11-2500

Digital ultrasonic cleaner is effective way to clean lab tools. It is a high efficiency cleaner specifically developed to clean surgical instruments.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

ULC11-2500 DIGITAL DESKTOP ULTRASONIC CLEANER

High strength and environmental ABS injection molding

Top quality ultrasonic ceramics parts, quality higher than competitors

Industrial IC, high voltage isolation design, safe using and long working life

Damp proof and anti corrosion process of the PCB

Strict production and quality test control

Digital control, multi stage ultrasonic cleaning adjustment

The tank is polished by SUS304 stainless steel stamping point

2500V high voltage test to ensure the product safety, strict production process, aging process and quality control

SPECIFICATIONS

Model	ULC11-2500
Capacity	2500 ml
Working frequency	40 Khz
Ultrasonic Power	70 W
Heating power	100 W
Timing time	1-8 mins
Temperature range	65°C
Cover	Yes
Drain valve	No
Washing basket	Standard
Inner tank dimension	247x150x78 mm
Packing size	325x245x210 mm
Gross weight	3.04 Kg
Power Supply	220V+ 10% 50Hz/60Hz 110V+ 10% 50Hz/60Hz

The logo for Labstac, featuring the word "LABSTAC" in a bold, sans-serif font. The letters "LAB" are in red and "STAC" are in black.

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com