

# DESKTOP MULTI FREQUENCY TYPE ULTRASONIC CLEANER ULC12-22



# DESKTOP MULTI FREQUENCY TYPE ULTRASONIC CLEANER ULC12-22

Desktop Multi frequency type Ultrasonic Cleaner uses multi frequency of 40/80/120 Khz for carrying out vigorous cleaning to remove dirt and heavy contaminants.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## ULC12-22 DESKTOP MULTI FREQUENCY TYPE ULTRASONIC CLEANER

Digital display and control

Special multi frequency ultrasonic transducer, cleaning conversion efficiency is higher

No stud welding technology, cleaning effect is better

MCU SWEEP multi frequency ultrasonic generator output, frequency sweep drive ultrasonic cleaning effect more powerful

MHC heating system, maximum temperature is up to 80°C

Whole machine noise reduction processing

All stainless steel material inner tank is made by punching SUS304 1.0mm

Multi frequency or three frequency cleaning modes can be automatically/ manually switched

2500V high voltage test to ensure the product safety, strict production, aging production quality control



### SPECIFICATIONS

Model	ULC12-22
Capacity	22 L
Working frequency	28/40 Khz or 40/80 Khz or 40/80/120 Khz
Ultrasonic Power	480 W
Heating power	500 W
Timing time	1-30 mins
Temperature range	0- 80°C
Cover	Yes
Drain valve	Yes
Washing basket	Optional
Inner tank dimension	500x300x150 mm
Packing size	657x454x430 mm
Gross weight	15 Kg
Power Supply	220V+ 10% 50Hz/60Hz 110V+ 10% 50Hz/60Hz



**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)