

# MECHANICAL ULTRASONIC DESKTOP CLEANER ULC14-11



# MECHANICAL ULTRASONIC DESKTOP CLEANER

## ULC14-11

Mechanical Ultrasonic desktop cleaner has high capacity for internal operations, high cleaning accuracy Effectively cleans small dirt particles, even deep holes, small seams and concealment of workpieces.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## ULC14-11 MECHANICAL ULTRASONIC DESKTOP CLEANER

Mechanical temperature and time setting

Industrial high Q transducers, cleaning conversion efficiency is higher

No stud- welding technology, Cleaning effect is better

MCU SWEEP multi frequency ultrasonic generator output, frequency sweep drive ultrasonic cleaning effect more powerful

MHC heating system, maximum temperature is up to 80°C

All stainless steel material inner tank is made by punching SUS304 1.0mm

2500V high voltage test to ensure the product safety, strict production, aging production quality control



### SPECIFICATIONS

Model	ULC14-11
Capacity	11 L
Working frequency	40 KHz
Ultrasonic Power	240 W
Heating power	300 W
Timing time	1- 60 mins
Temperature range	0- 80°C
Cover	Yes
Drain valve	Yes
Washing basket	Optional
Inner tank dimension	300x240x150 mm
Packing size	412x357x380 mm
Gross weight	8.5 kg
Power Supply	220V+ 10% 50Hz/60Hz 110V+ 10% 50Hz/60Hz

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)