

# MECHANICAL ULTRASONIC DESKTOP CLEANER ULC14-30A



# MECHANICAL ULTRASONIC DESKTOP CLEANER

## ULC14-30A

Mechanical Ultrasonic desktop cleaner has high capacity for internal operations, high cleaning accuracy Effectively cleans small dirt particles, even deep holes, small seams and concealment of workpieces.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## ULC14-30A MECHANICAL ULTRASONIC DESKTOP CLEANER

Mechanical temperature and time setting

Industrial high Q transducers, cleaning conversion efficiency is higher

No stud- welding technology, Cleaning effect is better

MCU SWEEP multi frequency ultrasonic generator output, frequency sweep drive ultrasonic cleaning effect more powerful

MHC heating system, maximum temperature is up to 80°C

All stainless steel material inner tank is made by punching SUS304 1.0mm

2500V high voltage test to ensure the product safety, strict production, aging production quality control



### SPECIFICATIONS

Model	ULC14-30A
Capacity	30 L
Working frequency	40 KHz
Ultrasonic Power	600 W
Heating power	600 W/220V 500 W/110V
Timing time	1- 60 mins
Temperature range	0- 80°C
Cover	Yes
Drain valve	Yes
Washing basket	Optional
Inner tank dimension	500x300x200 mm
Packing size	657x450x475 mm
Gross weight	16.0 Kg
Power Supply	220V+ 10% 50Hz/60Hz 110V+ 10% 50Hz/60Hz

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)