

MECHANICAL ULTRASONIC DESKTOP CLEANER ULC14-30B



MECHANICAL ULTRASONIC DESKTOP CLEANER

ULC14-30B

Mechanical Ultrasonic desktop cleaner has high capacity for internal operations, high cleaning accuracy Effectively cleans small dirt particles, even deep holes, small seams and concealment of workpieces.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

ULC14-30B MECHANICAL ULTRASONIC DESKTOP CLEANER

Mechanical temperature and time setting

Industrial high Q transducers, cleaning conversion efficiency is higher

No stud- welding technology, Cleaning effect is better

MCU SWEEP multi frequency ultrasonic generator output, frequency sweep drive ultrasonic cleaning effect more powerful

MHC heating system, maximum temperature is up to 80°C

All stainless steel material inner tank is made by punching SUS304 1.0mm

2500V high voltage test to ensure the product safety, strict production, aging production quality control



SPECIFICATIONS

| | |
|----------------------|---|
| Model | ULC14-30B |
| Capacity | 30 L |
| Working frequency | 28 Khz |
| Ultrasonic Power | 600 W |
| Heating power | 600 W/220V 500 W/110V |
| Timing time | 1- 60 mins |
| Temperature range | 0- 80°C |
| Cover | Yes |
| Drain valve | Yes |
| Washing basket | Optional |
| Inner tank dimension | 490x290x195 mm |
| Packing size | 657x450x475 mm |
| Gross weight | 20.7 Kg |
| Power Supply | 220V+ 10% 50Hz/60Hz 110V+ 10% 50Hz/60Hz |

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com